# LOXO data sheet LOXO data sheet LOXO data sheet LOXO

## Anti-C<sup>w</sup> (Anti-RH8) (Monoklonal)

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## Anti-C<sup>w</sup> (Anti-RH8) (Monoclonal)

For Slide and Tube Test  
FOR IN VITRO DIAGNOSTIC USE ONLY

### INTENDED USE

Monoclonal agglutinating Anti- C<sup>w</sup> - testserum is produced from cell culture supernatants of a heterohybridoma-cell line. The cells are secreting an antibody of IgM-type, that reacts specific with the corresponding antigen. The antibody is human protein. The testserum is used to determine whether red blood cells possess or lack the corresponding blood group antigen. The testserum is intended to be used by qualified and technical personnel only.

### PRINCIPLE OF PROCEDURE

The procedures used with this reagent is based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed toward the antigen.

### REAGENT

The listed reagent is produced by following cell clone:  
Anti-C<sup>w</sup> agglutinating (monoclonal, human IgM, clone MS-110)

This reagent contains <0.1% (w/v) sodium azide as preservative. Additionally the reagent is prepared of active antibody, sodium chloride, macromolecules and bovine albumin.

### WARNING

This reagents is prepared from supernatant of cell cultures. As biological products it should be looked upon as potentially infectious because of never complete exclusion of danger through excitants of disease. The reagent contains sodium azide, that may be toxic and may react with lead or copper to form highly explosive salts. Because of these reasons reagent should be handled with proper care.

### STORAGE REQUIREMENT

Store at 2 to 8°C. May be at room temperature (15 to 30°C) while in use. In principle, store and use the reagents to declared expiry date only.

### REMARKS

- Strength of positive reactions also depends on age of used blood
- With each testing positive and negative controls should be performed.
- Inappropriate storage impairs efficacy of the reagent.
- Centrifugation highly different from appointed relative centrifugal force may lead to false results
- Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at 2 to 8°C. Blood drawn into sodium citrate or EDTA should be tested within 14 days. Blood obtained by finger puncture may be tested directly by the slide method but, to avoid clotting, blood collected in this manner should be mixed quickly with the reagent.
- The procedures identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on automated systems.
- For usage of this testsera all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.

### REAGENT PREPARATION

There is no preparation of the reagent required. Use reagent directly from the vials.

### PROCEDURE

Not provided material, additionally needed

#### at Slide Method

- Glass slide
- Pasteur pipette
- Mixing stick

#### at Tube Centrifugation Method

- Test tubes, 10 x 75 mm or 12 x 75 mm
- Pipettes designed to deliver approximately 100 µL
- Centrifuge
- Isotonic saline (with 0,85 - 0,9% sodium chloride)

### Test procedure

#### Slide Method

- Use erythrocyte sediment or whole blood only.
- Place one drop (approximately 50 µL) of appropriate reagent on a glass slide.
- Using a Pasteur pipette add one drop erythrocyte sediment or whole blood (approximately 50 µL) to the glass slide.
- Mix the erythrocytes with reagent well with a stick and spread to a circle (diameter 2 cm)
- By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds). Unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.

#### Tube Centrifugation Method

- Use a 2% to 5% suspension of red blood cells in isotonic saline (cells washed one time or up to three times with isotonic saline) only.
- Add 100 µL (alternative: one drop = approximately 50 µL) of appropriate reagent to each tube
- Add 100 µL (alternative: one drop = approximately 50 µL) of appropriate cell suspension to each tube
- Mix well by slightly shaking.
- Incubate tube at room temperature (15 to 30 °C) for 15 min.
- Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 x g).
- Gently resuspend the red cells and check macroscopically for agglutination within 3 minutes. Document the result.

### INTERPRETATION OF RESULTS

" Slightly rotating / shaking " at Slide Method and at Tube Centrifugation Method  
Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.  
Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

### LIMITATIONS OF THE PROCEDURE

- Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
- No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
- Enzyme treated erythrocytes may react unspecific.
- Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may give weaker reactivity compared to control cells.
- Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test (DAT)) may give weak reactions. In extreme cases, false-negative results may occur.

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