LOXO data sheet LOXO data sheet LOXO data sheet LOXO

Anti-D blend (monoklonal, human IgG + IgM)

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FOR IN VITRO DIAGNOSTIC USE ONLY

INTENDED USE

Monoclonal Anti-D blend - testserum is produced from cell culture supernatants of hetero-hybridoma cell lines. One cell line is secreting an antibody of IgG-type and the other an antibody of IgM-type, that both react specific with the corresponding antigen. The testserum is used to determine whether human red blood cells possess or lack the corresponding blood group antigen D. This testserum reacts negatively with erythrocytes of category D⁺ with the IgM-antibody and positive with the IgG-antibody in indirect Coombs-test. The testserum is intended to be used by qualified and technical personnel only.

PRINCIPLE OF PROCEDURE

The procedures used with this reagent are based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed towards the antigen.

REAGENTS

The listed reagent contains antibodies of the following cell clones:

Anti-D blend (monoclonal, human IgG + IgM, clones: MS-26; TH-28)

The reagent contains <0.1% (w/v) sodium azide as preservative. Additionally the reagent is prepared of active antibody, sodium chloride, macromolecules and bovine albumin.

WARNING

This reagent was prepared from supernatants of cell cultures. As a biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through extants of disease. The reagent contains sodium azide, that may be toxic and may react with lead or copper to form highly explosive salts. Because of these reasons this reagent should be handled with proper care.

STORAGE REQUIREMENT

Store at 2 to 8 °C. May be at room temperature (15 to 30 °C) while in use. In principle, store and use the reagents to indicated expiry date only.

REMARKS

- Strength of positive reactions also depends on age of used blood
- With each testing positive and negative controls should be performed.
- Inappropriate storage impairs efficacy of the reagent.
- Centrifugation highly different from appointed relative centrifugal force may lead to false results.
- Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at 2 to 8 °C. Blood drawn into sodium citrate or EDTA should be tested within 14 days. Blood obtained by finger puncture may be tested directly by the slide method but, to avoid clotting, blood collected in this manner should be mixed quickly with the reagent.
- The procedures identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on automated systems.
- For usage of this testsera all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.

REAGENT PREPARATION

There is no preparation of the reagents required. Use reagents directly from the vials.

PROCEDURE

Not provided material, additionally needed

at Slide Method:

Glass slide; Pasteur pipette; Mixing stick

at Tube Centrifugation Method:

Test tubes, 10 x 75 mm or 12 x 75 mm; Pipettes designed to deliver approximately 100 µL; Centrifuge; Isotonic saline (with 0,85 - 0,9% sodium chloride)

Test procedure

Slide Method

- Use erythrocyte sediment or whole blood only.
- Place one drop (approximately 50 µL) of appropriate reagent on a glass slide.
- Using a Pasteur pipette add one drop erythrocyte sediment or whole blood (approximately 50 µL) to the glass slide.
- Mix the erythrocytes with reagent well with a stick and spread to a circle (diameter 2 cm)
- By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds). Unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.

Tube Centrifugation Method

- Use a 2% to 5% suspension of red blood cells in isotonic saline (cells washed one time or up to three times with isotonic saline) only.
- Add 100 µL (alternative: one drop = approximately 50 µL) of appropriate reagent to each tube.
- Add 100 µL (alternative: one drop = approximately 50 µL) of appropriate cell suspension to each tube.
- Mix well by slightly shaking.
- Incubate tube at room temperature (15 to 30 °C) for 1 - 15 min.
- Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800-1.000 x g).
- Enzyme treated red cells and check macroscopically for agglutination within 3 minutes. Document the result.

Indirect Coombs-test

- Repeat step 1.-4. with fresh material or continue after reading of result directly.
- Incubate tube at 37 °C for 30 min.
- Wash red cells 3 times with (cold) isotonic saline.
- Add 100 µL Anti-Human-Globulin-serum (Coombs-serum) to each tube
- Centrifugation of tube for 1 min at 1.000 rpm (approximately 180-270 x g).
- Gently resuspend the red cells and check macroscopically for agglutination within 3 minutes. Document the result.

INTERPRETATION OF RESULTS

" Slightly rotating / shaking " at Slide Method / at Tube Centrifugation Method

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

LIMITATIONS OF THE PROCEDURE

- Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
- No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
- Enzyme treated erythrocytes may react unspecific.
- Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may give weaker reactivity compared to control cells.
- Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test (DAT)) may give weak reactions. In extreme cases, false-negative results may occur.
- With slide method most of weak D-antigens (Dweak) and categories are not recognized.
- This testserum reacts negatively with erythrocytes of category DVI in an agglutinating test and positively in an indirect Coombs-test.
- The described procedures for testing are for usage with Coombs-serum produced by LOXO GmbH. On principle Coombs-sera of other manufacturers can be used, but laboratories have to follow the procedures of this manufacturer. Also laboratories must follow approved validation procedures to demonstrate compatibility of this product.
- If products are taken out of 1 L flasks, the material of this bottle has to be used within 6 month.

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