# LOXO data sheet LOXO data sheet LOXO data sheet LOXO

## Anti-Le<sup>a</sup> (Anti-LE1) (Monoklonal) Anti-Le<sup>b</sup> (Anti-LE2) (Monoklonal)

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NUR ZUR IN-VITRO-DIAGNOSTIC

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FOR IN VITRO DIAGNOSTIC USE ONLY

### INTENDED USE

Monoclonal agglutinating Anti-Le<sup>a</sup> and Anti-Le<sup>b</sup>-testsera are produced from cell culture supernatants of hybridoma-cell lines. The cells are secreting antibodies of IgM-type, that reacts specific with the corresponding antigen. The testsera are used to determine whether red blood cells possess or lack the corresponding blood group antigen. The testsera are intended to be used by qualified and technical personnel only.

### PRINCIPLE OF PROCEDURE

The procedures used with this reagent are based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed toward the antigen.

### REAGENTS

The listed reagent is produced from following cell clone:

- Anti-Le<sup>a</sup> agglutinating (monoclonal, clone 78FR2.3)
- Anti-Le<sup>a</sup> agglutinating (monoclonal, clones LM 112/161)
- Anti-Le<sup>b</sup> agglutinating (monoclonal, clones LEA2)
- Anti-Le<sup>b</sup> agglutinating (monoclonal, clone 96FR2.10)
- Anti-Le<sup>a</sup> agglutinating (monoclonal, Klon LM129/181)
- Anti-Le<sup>b</sup> agglutinating (monoclonal, Klon LEB2)

The reagent contains <0.1% (w/v) sodium azide as preservative. Additionally the reagent is prepared of active antibody, sodium chloride, macromolecules and bovine albumin.

### WARNING

These reagents were prepared from supernatants of cell cultures. As biological products it should be looked upon as potentially infectious because of never complete exclusion of danger through excipients of disease. These reagents contain sodium azide, that may be toxic and may react with lead or copper to form highly explosive salts. Because of these reasons the reagents should be handled with proper care.

### STORAGE REQUIREMENT

Store at 2 to 8°C. May be at room temperature (15 to 30°C) while in use. In principle, store and use the reagents to declared expiry date only.

### REMARKS

- Strength of positive reactions also depends on age of used blood
- With each testing positive and negative controls should be performed.
- Inappropriate storage impairs efficacy of the reagent.
- Centrifugation highly different from appointed relative centrifugal force may lead to false results
- Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at 2 to 8°C. Blood drawn into sodium citrate or EDTA should be tested within 14 days.
- The procedures identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on automated systems.
- For usage of this testsera all effective national laws, directives and guidelines have to be observed, in Germany especially the „Richtlinien zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Hämotherapie)“.

### REAGENT PREPARATION

There is no preparation of the reagents required. Use reagents directly from the vial.

### PROCEDURE

Not provided material, additionally needed

#### at Tube Centrifugation Method

- Test tubes, 10 x 75 mm or 12 x 75 mm
- Pipettes designed to deliver approximately 100 µL
- Centrifuge
- Isotonic saline (with 0,85 - 0,9% sodium chloride)

#### Test procedure

#### Tube Centrifugation Method

- Use a 2% to 5% suspension of red blood cells in isotonic saline (cells washed one time or up to three times with isotonic saline) only.
- Add 100 µL (alternative: one drop = approximately 50 µL) of appropriate reagent to each tube
- Add 100 µL (alternative: one drop = approximately 50 µL) of appropriate cell suspension to each tube
- Mix well by slightly shaking.
- Incubate tube at room temperature (15 to 30 °C) for 15 min.
- Centrifugation of tube for 1 minute at 1.000 rpm (approximately 180-270 x g).
- Gently resuspend the red cells and check macroscopically for agglutination within 3 minutes. Document the result.

#### INTERPRETATION OF RESULTS

"Slightly shaking" at Tube Centrifugation Method

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

#### LIMITATIONS OF THE PROCEDURE

- Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
- No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
- Enzyme treated erythrocytes may react unspecific.
- Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may give weaker reactivity compared to control cells.
- Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test (DAT)) may give weak reactions. In extreme cases, false-negative results may occur.

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