

# Pefachrome<sup>®</sup> PL-Strept

**Application:** Highly sensitive chromogenic peptide substrate for the plasmin-streptokinase complex. Determination of plasminogen levels for research, in-process and quality control.

**Formula:** H-D-Nle-CHA-Arg-pNA·2AcOH

**Principle:** Plasminogen + streptokinase ==> plasminogen-streptokinase complex (PSC)  
H-D-Nle-CHA-Arg-pNA + PSC ==> H-D-Nle-CHA-Arg-OH+ pNA + PSC

**K<sub>M</sub>:** 0.4 mM **V<sub>max</sub>:** 0.024 μmol/min

**MW:** 680.8

**Storage:** May be used by the expiry date given on the label when stored unopened, protected from moisture, in the dark, 2-8°C. Avoid contamination of the reagents by micro-organisms.  
Shipment of product does not require cooling during the time of transportation.

**Material required but not provided:**

Buffer, reference material

**Buffer:** 50 mM Tris-imidazole buffer pH 7.2, 200 mM NaCl

**Assay:** Suggested protocol for the determination of plasminogen-streptokinase-complex:

0.100 ml citrated plasma, diluted 1 : 10 in buffer  
0.100 ml streptokinase, 10'000 IU/ml in dist. water  
=> incubate for 3 min at 37° C  
1.600 ml buffer  
0.200 ml Pefachrome<sup>®</sup> PL-Strept (2 mM in water)  
=> Determination of ΔOD/min at 405 nm

**Reference:** Svendsen LG, Fareed J, Walenga JM, Hoppensteadt D.  
Newer synthetic peptide substrates in coagulation testing: Some practical considerations for automated Methods. Semin Thromb Haemost 1983; 9: 250-62.

**Package size:** Vial containing 25 mg  
Bulk [g]

**Code:** 083-05

FOR RESEARCH USE ONLY. NOT FOR HUMAN USE OR DRUG USE.