# **IMMBIOMED**

#### ImmBioMed GmbH & Co. KG

Bergstrasse 85, DE-64319 Pfungstadt Phone: +49 6157 9799 685, Fax: +49 6157 86443

# Synthetic Inhibitor of Fibrin Polymerization

Product No. ADG4530N

## **Description**

The synthetic inhibitor binds with a high affinity to fibrinogen, inhibits fibrin polymerization, modifies the mechanical properties of fibrin clots and can dissociate non-stabilized fibrin gels.

## Chemistry

Formula: H-Gly-Pro-Arg-Pro-OH×AcOH

**MW:** 485.5

Solubility: readily soluble in water

## **Application**

The inhibitor is added to fibrinogen containing reaction mixtures to inhibit disturbing fibrin-related turbidity, gel formation and fibrin deposition in diagnostic and preparative procedures (e.g. measurement of thrombin generation in plasma). The inhibitor is also used to inhibit fibrin formation during purification and processing of clotting factors and other plasma proteins.

#### Recommended inhibitor concentrations:

Application	Concentration
Gelation prevention:	
des-AB-fibrin	0.5-2.3 mM
des-A-fibrin	2.3 mM
human plasma, HC-II test	0.5 mg/ml
bov. fibrinogen 0.2% / thrombin 0.5 U/ml	1.5 mg/ml
Clot liquefaction:	
des-AB-fibrin	15.0 mM
des-A-fibrin	1.0 mM

#### **Presentation**

Vial containing 50 mg lyophilized inhibitor.

#### Reconstitution

Reconstitute with filtered deionized water.

# Storage and Stability

May be used by the expiry date given on the label when stored unopened, protected from moisture, in the dark, 2-8°C. Avoid contamination of the reagents by micro-organisms. Shipment of product does not require cooling during the time of transportation.

#### References

Effect of calcium and synthetic peptides on fibrin polymerization. Furlan M et al., Thromb Haemost 1982; 47: 118-21.

The ability of thrombin inhibitors to reduce the thrombin activity generated in plasma on extrinsic and intrinsic Activation. Prasa Det al., Thromb Haemost 1997: 77: 498-503.



IMMUNOLOGIE • MOLEKULARBIOLOGIE BIOCHEMIE • PRODUKTE UND SYSTEME

Gerhart-Hauptmann-Str. 48 69221 Dossenheim

Tel +49 6221 868023 Fax +49 6221 8680255 www.loxo.de - info@loxo.de

#### Hinweis/Note:

Der Packungsbeileger dient nur als erste Information. Der relevante Packungsbeileger liegt der Ware bei.

The datasheet is for information purposes only. The current datasheet will be enclosed with product shipment.