IMMBIOMED

ImmBioMed GmbH & Co. KG

Bergstrasse 85, DE-64319 Pfungstadt Phone: +49 6157 9799 685. Fax: +49 6157 86443

Factor X Deficient Plasma

Product No. ADG210

Description

Factor X Deficient Plasma is manufactured from pooled normal human plasma depleted of Factor X using antibodies directed to Factor X. Only the highest quality citrated plasma is used as starting material.

Specification

The plasma is assayed to ensure the proper activity of the remaining factors (>50% activity), and is tested for fibrinogen, PT, aPTT. *Please note: Results may vary depending on materials and methods used.

Activity assay	Lot xxxxxx *
PT	xxx sec
aPTT	xxx sec
Fibrinogen	xxx mg/dL
Factor II	xx %
Factor V	xx %
Factor VII	xx %
Factor VIII	xx %
Factor IX	xx %
Factor X	< 1 %
Factor XI	xx %
Factor XII	xx %

Application

For research assays of Factor X antigen and functionality as a negative control and for dilution of standards and test samples.

Warning

The plasma is prepared from screened normal citrated human plasma that have been tested and found negative for HIV and Hepatitis by FDA-approved methods.

As no known test method can provide complete assurance that products derived from human blood will not transmit HBsAg, HCV, HIV-1, HIV-2 or other blood-borne pathogens, this reagent should be handled as recommended for any potentially infectious human specimen.

Presentation

Set of 10 clear, screw-capped vials each containing 1 mL of frozen plasma.

Preparation and Stability

Thaw 1 ml vials in 37 °C water bath for 5 minutes. Product is stable until date stated on vial label when stored at -20 °C. Once thawed, plasma is stable for 4 hours at 2-8 °C in original vial.

Storage

Plasma is shipped frozen and should be stored at -20 °C or colder.

Distributed by:

LOXO

IMMUNOLOGIE • MOLEKULARBIOLOGIE BIOCHEMIE • PRODUKTE UND SYSTEME

Gerhart-Hauptmann-Str. 48 69221 Dossenheim

Tel +49 6221 868023 Fax +49 6221 8680255 www.loxo.de - info@loxo.de

Hinweis/Note:

Der Packungsbeileger dient nur als erste Information. Der relevante Packungsbeileger liegt der Ware bei.

The datasheet is for information purposes only. The current datasheet will be enclosed with product shipment.