IMMBIOMED

ImmBioMed GmbH & Co. KG

Bergstrasse 85, DE-64319 Pfungstadt Phone: +49 6157 9799 685, Fax: +49 6157 86443

Monoclonal antibodies against Streptococcal proteinase Product Nos.

ADG0303L, ADG0304L, ADG0305L, ADG0306L and ADG0310L

Description

Streptococcal proteinase is a sulfhydryl enzyme elaborated by group A streptococci. It is excreted into the medium as a zymogen (44 kDa). The inactive enzyme (32 kDa) formed by proteolysis of the zymogen requires activation by exposure to reducing agents.

The reduced, active enzyme can be readily inactivated by reagents known to react with sulfhydryl groups.

Preparation

The monoclonal antibodies were generated against purified Streptococcal proteinase. The antibodies have been purified from cell culture supernatant using Protein G affinity chromatography.

Presentation

Screw capped vial containing 1 mg of purified antibody in PBS pH 7.4. The IgG concentration is given on the vial label. Spin the vial briefly before opening.

Storage and Stability

Store the antibody at $2^{\circ}-8^{\circ}$ C. For long-term storage the antibody should be aliquoted and stored at -20° C or colder. It is recommended to avoid freeze-thaw cycles.

Known applications

Western blot

The antibodies are suitable for western blot analysis. They show different reactivity towards the zymogen (Z) or the active proteinase (ST). The clones HD-SP 3.1, 4.1 and 5.1 react solely with the zymogen. Apparently, they recognize an epitope on the 13K-fragment, which is cleaved of during proteolytic activation.

Product no.	clone	Z	ST
ADG0303L	HD-SP 3.1	+	-
ADG0304L	HD-SP 4.1	+	-
ADG0305L	HD-SP 5.1	+	-
ADG0306L	HD-SP 6.1	+	+
ADG0310L	HD-SP 10.1	-	+

Distributed by:



IMMUNOLOGIE • MOLEKULARBIOLOGIE BIOCHEMIE • PRODUKTE UND SYSTEME

Gerhart-Hauptmann-Str. 48 69221 Dossenheim

Tel +49 6221 868023 Fax +49 6221 8680255 www.loxo.de - info@loxo.de

Hinweis/Note:

Der Packungsbeileger dient nur als erste Information. Der relevante Packungsbeileger liegt der Ware bei. The datasheet is for information purposes only. The current datasheet will be enclosed with product shipment.