

## Fibrinogen Deficient Plasma

Product No. ADG204

Lot XXXXXX

### Description

Fibrinogen Deficient Plasma is manufactured from pooled normal human plasma defibrinated under controlled conditions, using purified human thrombin. Only the highest quality citrated plasma is used as starting material.

### Specification

Activity assay	Lot XXXXXX
PT (Stago Neoplastine C1+)	>XXX sec
APTT (HemosIL APTT-SP)	>XXX sec
Factor II (activity)	XXX U/ml
Factor V (activity)	XXX U/ml
Factor VIII (activity)	XXX U/ml
Fibrinogen (antigen) by ELISA	XXX g/L

### Application

For research assays of fibrinogen antigen and functionality, as a negative control and for dilution of standards and test samples.

### Distributed by:



IMMUNOLOGIE • MOLEKULARBIOLOGIE  
BIOCHEMIE • PRODUKTE UND SYSTEME

Gerhart-Hauptmann-Str. 48  
69221 Dossenheim

Tel +49 6221 868023  
Fax +49 6221 8680255

www.loxo.de - info@loxo.de

### Warning

The plasma is prepared from normal citrated human plasma collected from donors screened for CJD and which was tested and found negative for HBsAg, syphilis (RPR or T Pallidum-G-EIA) and antibodies to HIV and HCV and non-reactive for HIV-1 RNA and HCV RNA by FDA approved methods.

As no known test method can provide complete assurance that products derived from human blood will not transmit HBsAg, HCV, HIV-1, HIV-2 or other blood-borne pathogens, this reagent should be handled as recommended for any potentially infectious human specimen.

### Presentation

Set of 5 clear, screw-capped vials each containing 1 mL frozen plasma. The plasma contains 20mM HEPES.

### Preparation and Stability

Thaw 1 ml vials in 37 °C water bath for 5 minutes. Product is stable until date stated on vial label when stored at -60 °C. Once thawed, plasma is stable for 4 hours at 2-8 °C in original vial.

### Storage

Plasma is shipped frozen and should be stored at -60 °C.

### Hinweis/Note:

Der Packungsbeileger dient nur als erste Information. Der relevante Packungsbeileger liegt der Ware bei.  
The datasheet is for information purposes only. The current datasheet will be enclosed with product shipment.

For research use only!

ADG204©08072019\_Lot