# **IMMBIOMED**

#### ImmBioMed GmbH & Co. KG

Bergstrasse 85, DE-64319 Pfungstadt Phone: +49 6157 9799 685, Fax: +49 6157 86443

# Monoclonal antibody against Outer surface protein A (OspA) Borrelia burgdorferi [LA84.1]

Product No. ADG0080L

# **Description**

The outer surface Protein A (OspA) of *Borrelia burgdorferi* has been a focus of scientific scrutiny since its identification. This 31-kDa lipoprotein is localized on the surface and in the periplasmic space of *B. burgdorferi*. OspA coats the spirochetes from the time the enter ticks until they are transmitted to a mammal. Vaccination with OspA elicits antibody that can target spirochetes in the tick midgut during feeding and inhibit transmission to mammals. OspA was the primary component of the human LYMErix™ vaccine.

## **Properties**

The monoclonal antibody ADG0080L (**clone LA84.1**) is a murine monoclonal antibody, subclass IgG<sub>2a</sub> recognizing OspA. Mice were immunized with cell lysates of *Borrelia burgdorferi*. The antibody has been purified from cell culture supernatant using Protein G affinity chromatography.

#### **Presentation**

Vial containing 1 mg purified antibody in PBS pH 7.4. The concentration is given on the vial label. Spin the vial briefly before opening.

# Storage and Stability

Store the antibody at 2°-8°C. For long-term storage the antibody should be aliquoted and stored at -20°C or colder. It is recommended to avoid freeze-thaw cycles.

### **Applications**

#### A. ELISA

The antibody can be used as capture antibody in ELISAs. An antibody concentration of 1-10 µg/ml is recommended.

# **B.** Westernblot

The antibody is suitable for Western blot analysis, detecting native and recombinant OspA following SDS-PAGE under reducing conditions. A primary antibody concentration of 1-10 µg/mL is recommended.

# C. Immunocytochemistry

The antibody can be used for immunocytochemistry on paraformaldehyde fixed spirochetes and for cryoelectron tomography.

#### References

- Cloning and sequencing of the gene endcoding the outer surface prtotein A (OspA) of a European Borrelia burgdorferi isolate. Wallich et al. Nucleic Acids Res. 1989; 17(21):8864
- Monoclonal antibodies specific for the Outer surface protein A (OspA) of *Borrelia burgdorferi* prevent Lyme borreliosis in severe combined immunodeficiency (scid) mice. Schaible et al. *Proc. Natl. Acad. Sci. USA* 1990; 87(10):3768-3772
- Vaccination against Lyme disease with recombinant Borrelia burgdorferi outer-surface lipoprotein A with adjuvant. Lyme Disease Vaccine Study Group. Steere et al. N. Engl. J. Med. 1998; 339(4):209-215
- Structural identification of a key protective B-cell epitope in Lyme disease antigen OspA. Ding et al. *J. Mol. Biol.* 2000; 302(5):1153-1164
- 5. Nanoscopic localization of surface-exposed antigens of *Borrelia burgdorferi*. Lemgruber et al. *Microsc. Microanal*. 2015; 21(03):680-688
- 6. Comparative cryo-electron tomography of pathogenic Lyme disease spirochetes. Kudryashev et al. Mol. Microbiol. 2009; 71(6):1415-1434



Tel +49 6221 868023 Fax +49 6221 8680255 www.loxo.de - info@loxo.de

For research use only!